

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

Bird of Prey genes in the form of CloneCards or CloneTabs can be stored at room temperature (15-25°C) for up to 18 months.

For more information, please refer to the Bird of Prey Handbook, which can be found at www.ospreybio.com. For technical assistance, please contact us at clonecard@ospreybio.com

Notes before starting

- It is recommended to perform a plasmid prep for each gene involved in your vector design to ensure enough DNA is available for subcloning.
- We recommend 1µl of enzyme in restriction digests regardless of units/µl to ensure consistency, efficiency, and proper enzyme activity.
- Never use the same Transcription Unit (TU) or IRES element in a single vector construct. This will lead to vector recombination.
- If a 3rd Effector is desired, please refer to the Optional section beginning on Page 4. Note that you must still start at Step 1 below.

Materials Needed:

- BoP Vectors: Effectors, IRES elements, Transcription Units
- Restriction enzymes: EcoRV, PacI, PvuI, KsaI, Scal, XhoI, NotI, Sall
- T4 DNA Ligase and buffer
- Gel electrophoresis equipment
- DNA purification kit
- Competent E. coli cells

Step-by-Step Procedure:

1. Design

- Create a table with chosen design combinations

2. Digest I – First Effector

- Set up two separate digestion reactions:
 - a) Donor Vector (IRES): PacI + KsaI
 - b) Insert Vector (Effector/ORF): PvuI + KsaI
- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Sample ID	Volume (µL)							CIP	H ₂ O
	Buffer	Vector Acceptor (pDNA Backbone)	Insert Donor (DNA)	Enzyme I: PacI	Enzyme II: KsaI	Enzyme III: PvuI			
IRES-1	4	20	-	1	1	-	1	13	
IRES-2	4	20	-	1	1	-	1	13	
E1	4	20	-	-	1	1	-	13	
E3	4	20	-	-	1	1	-	13	

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

3. Gel Electrophoresis

- Run digested samples on a 1% agarose gel using the above table volumes
- Visualize bands under UV light

4. Fragment Purification

- Excise desired bands from gel (try using [OspreyBio's gel band cutters](#) for faster and more accurate results)
- Purify DNA using gel extraction kit (we recommend Zymo's "Zymoclean Kit" – Cat #D400 series)

5. Ligation

- Mix purified Controller and Insert fragments in individual and labeled tubes
- Add T4 DNA Ligase and buffer (see below table for values)
- Incubate at room temperature overnight.
- After 12-24 hours, heat shock at 37°C to stop reaction

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor (pDNA Backbone)	Insert	Ligase	H ₂ O
IRES-1 (control)	10	5	-	1	84
IRES-1 + E1	10	5	10	1	74
IRES-2 (control)	10	5	-	1	84
IRES-2 + E3	10	5	10	1	74

6. Digest II – Second Effector

Sample ID	Volume (µL)							
	Buffer	Vector Acceptor (pDNA Backbone)	Insert Donor (DNA)	Enzyme I: EcoRV	Enzyme II: PvuI	Enzyme III: ScaI	CIP	H ₂ O
IRES-1-E1	4	20	-	1	1	-	1	13
E2	4	-	20	-	1	1	-	14
IRES-2-E3	4	20	-	1	1	-	1	13
E4	4	-	20	-	1	1	-	14

Digest:

- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Gel Electrophoresis: Perform in the same manner as above.

Fragment Purification: Perform in the same manner as above.

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

7. Ligation II

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor	Insert	Ligase	H ₂ O
IRES-1-E1 (control)	10	5	-	1	84
E2 + IRES-1-E1	10	5	10	1	74
IRES-2-E3	10	5	-	1	84
E4 + IRES-2-E3	10	5	10	1	74

8. Transformation

- Transform ligation mixture into competent E. coli
- Plate on Kanamycin antibiotic-containing agar
- Incubate overnight at 37°C

9. Pick Colonies & perform miniprep

- Pick two colonies from incubated plates
- Add each to a culture tube of 3ml kanamycin LB broth
- Place in a shaker incubator at 37°C overnight
- Perform miniprep to isolate plasmid DNA
- Verify construct via DNA sequencing
- Proceed with downstream applications (e.g., protein expression)

10. Digest III (IV if 3-Effectors) – Multi-Effectors into Transcription Units

Sample ID	Volume (µL)						
	Buffer	Vector Acceptor	Insert Donor	Enzyme I: PvuI	Enzyme II: KsaI	CIP	H ₂ O
TU-1	4	20	-	1	1	1	13
E2-IRES-1-E1	4	-	20	1	1	0	14
TU-2	4	20	-	1	1	1	13
E4-IRES-2-E3	4	-	20	1	1	0	14

Digest:

- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Gel Electrophoresis: Perform in the same manner as above.

Fragment Purification: Perform in the same manner as above.

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

11. Ligation III (IV if 3-Effectors)

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor	Insert	Ligase	H ₂ O
TU-1 (control)	10	5	-	1	84
E2-IRES-1-E1 + TU-1	10	5	10	1	74
TU-2 (control)	10	5	-	1	84
E4-IRES-2-E3 + TU-2	10	5	10	1	74

12. Sequence Verification

Run transformation protocol, pick colonies, & perform minipreps to verify sequence (steps on Page 3).

13. Digest IV (V if 3-Effectors)

Sample ID	Volume (µL)							
	Buffer	Vector Acceptor	Insert Donor	Enzyme I: XhoI	Enzyme II: NotI	Enzyme III: Sall	CIP	H ₂ O
TU-1	4	20	-	1	1	-	1	13
TU-2	4	-	20	-	1	1	0	14

Digest:

- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Gel Electrophoresis: Perform in the same manner as above.

Fragment Purification: Perform in the same manner as above.

14. Ligation IV (V if 3-Effectors)

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor	Insert	Ligase	H ₂ O
TU-1 (control)	10	5	-	1	84
TU-1 + TU-2	10	5	10	1	74

15. Sequence Verification

Run transformation protocol, pick colonies, & perform minipreps to verify sequence (steps on Page 3).

Remember to always refer to the full Bird of Prey handbook for detailed protocols and troubleshooting advice. This QuickStart guide is intended as a concise reference for experienced users.

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

OPTIONAL: Three-Effector Vectors

If a third ORF is desired, a second IRES must be used. We recommend limiting each Transcription unit to three ORFs to mitigate against poor expression levels.

If after designing the construct you wish to create, a third ORF is needed, we recommend performing the restriction digest of the desired Effectors and second IRES at the same time as the first IRES and its first ORF to save time.

Digest I – Simultaneous IRES-1 & IRES-2

Please note that the enzymes will be the same for both sets of IRESes.

- a. IRESes: PacI (E₁) and KasI (E₂).
- b. Effectors/ORFs: PvuI (E₁) and KasI (E₂)

Sample ID	Volume (µL)						
	Buffer	Vector Acceptor	Enzyme I: PacI	Enzyme II: KasI	Enzyme III: PvuI	CIP	H ₂ O
IRES-1	4	20	1	1	-	1	13
E1	4	20	-	1	1	0	14
IRES-2	4	20	1	1	-	1	13
E3	4	20	-	1	1	0	14
IRES-3	4	20	1	1	-	1	13
E5	4	20	-	1	1	0	14
IRES-4	4	20	1	1	-	1	13
E6	4	20	-	1	1	0	14

Digest:

- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Gel Electrophoresis: Perform in the same manner as above.

Fragment Purification: Perform in the same manner as above.

Bird of Prey™ QUICKSTART: Multi-Effector Multigene

Ligation I

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor	Insert	Ligase	H ₂ O
IRES-1 (control)	10	5	-	1	84
IRES-1 + E1	10	5	10	1	74
IRES-2 (control)	10	5	-	1	84
IRES-2 + E3	10	5	10	1	74
IRES-3 (control)	10	5	-	1	84
IRES-3 + E5	10	5	10	1	74
IRES-4 (control)	10	5	-	1	84
IRES-4 + E6	10	5	10	1	74

Digest & Ligation II

Proceed with second Digest & Ligation for the first IRES to make your 2-Effector IRES (Steps 6 & 7).

Digest III

Sample ID	Volume (µL)							CIP	H ₂ O
	Buffer	Vector Acceptor	Insert Donor	Enzyme I: Scal	Enzyme II: KasI	Enzyme III: EcoRV			
E2-IRES-1-E1	4	20	-	1	1	-	1	13	
IRES-3-E5	4		20		1	1	0	14	
E4-IRES-2-E3	4	20	-	1	1	-	1	13	
IRES-4-E6	4		20		1	1	0	14	

Digest:

- Incubate at 37°C for 2 hours
- Add CIP and incubate at 37°C for 1 hour

Gel Electrophoresis: Perform in the same manner as above.

Fragment Purification: Perform in the same manner as above.



Bird of Prey™ QUICKSTART: Multi-Effector Multigene

Ligation III

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor	Insert	Ligase	H ₂ O
E2-IRES-1-E1 (control)	10	5	-	1	84
E1-IRES-1-E2+IRES-2-E5	10	5	10	1	74
E4-IRES-2-E3 (control)	10	5	-	1	84
E1-IRES-1-E2+IRES-4-E6	10	5	10	1	74

Sequence Verification

Run transformation protocol, pick colonies, & perform minipreps to verify sequence (steps on Page 3).

Digest IV & Ligation IV

Proceed to Step 10 to subclone your 3-Effector vector into your Transcription Unit of choice and proceed with the rest of the protocol.