

## Bird of Prey™ QUICKSTART: Mono-Effector Monogene

Bird of Prey genes in the form of CloneCards or CloneTabs can be stored at room temperature (15-25°C) for up to 24 months.

For more information, please refer to the Bird of Prey Handbook, which can be found at [www.ospreybio.com](http://www.ospreybio.com). For technical assistance, please contact us at [clonecard@ospreybio.com](mailto:clonecard@ospreybio.com)

### Notes before starting

- It is recommended to perform a plasmid prep for each gene involved in your vector design to ensure enough DNA is available for subcloning.
- We recommend 1µl of enzyme in restriction digests regardless of units/µl to ensure consistency, efficiency, and proper enzyme activity.

### Materials Needed:

- BoP Vectors: Effectors & Controllers
- Restriction enzymes: PvuI & KasI
- T4 DNA Ligase and buffer
- Quick CIP
- Gel electrophoresis equipment
- DNA fragment purification kit
- Competent E. coli cells

## Step-by-Step Procedure:

### 1. Design

- Choose desired BoP Transcription Unit(s) and Effector vector(s).

### 2. Digest

- Set up two separate digestion reactions:
  - a) Controller Vector: PvuI + KasI
  - b) Insert Vector (Effector/ORF): PvuI + KasI
- Incubate for 2 hours at 37°C
- Add CIP and incubate at 37°C for 1 hour

Sample ID	Volume (µL)						
	Buffer	Vector Acceptor (pDNA Backbone)	Insert Donor (DNA)	PvuI	KasI	CIP	H <sub>2</sub> O
TU-1	4	20	-	1	1	1	13
E1	4	-	20	1	1	-	14

### 3. Gel Electrophoresis

- Run digested samples on a 1% agarose gel using the above table volumes (in µl).
- Visualize bands under UV light.

### 4. Fragment Purification

- Excise desired bands from gel (try using [OspreyBio's gel band cutters](#) for faster and more accurate results).
- Purify DNA using gel extraction kit (we recommend Zymo's "Zymoclean Kit" – Cat #D400 series).

## Bird of Prey™ QUICKSTART: Mono-Effector Monogene

---

### 5. Ligation

- Mix purified Controller and Insert fragments in individual and labeled tubes.
- Add T4 DNA Ligase and buffer (see below table for values).
- Incubate at room temperature overnight.
- After 12-24 hours, heat shock at 37°C to stop reaction

Sample ID	Volume (µL)				
	Buffer	Vector Acceptor (pDNA Backbone)	Insert	Ligase	H <sub>2</sub> O
<b>TU-1 (control)</b>	10	5	-	1	84
<b>TU-1 + E1</b>	10	5	10	1	74

### 6. Transformation

- Transform ligation mixture into competent E. coli
- Plate on Kanamycin antibiotic-containing agar
- Incubate overnight at 37°C

### 7. Pick Colonies & perform miniprep

- Pick two colonies from incubated plates
- Add each to a culture tube of 3ml kanamycin LB broth
- Place in a shaker incubator at 37°C overnight
- Perform miniprep to isolate plasmid DNA
- Verify construct via DNA sequencing
- Proceed with downstream applications (e.g., protein expression)

Remember to always refer to the full Bird of Prey handbook for detailed protocols and troubleshooting advice. This QuickStart guide is intended as a concise reference for experienced users.